


Anti-Eif3A Polyclonal Antibody For Wb, If/Icc, Elisa - Q14152

Item Number: CM0018467



Introduction

Rabbit polyclonal antibody targeting human EIF3A, a subunit of the eIF3 complex crucial for translation initiation. Validated in WB, IF/ICC, ELISA; cross-reacts with mouse. Ideal for research on protein synthesis, cell proliferation, and viral infection.

[Learn More](#)

Parameter	Value
Product Name	EIF3A Rabbit pAb
Remarks/Alias	EIF3; P167; p180; p185; TIF32; EIF3S10; eIF3-p170; eIF3-theta; EIF3A
Species	Human
GeneID (Human)	8661
GeneID	8661
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human EIF3A (NP_003741.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14152
Protein Weight	167kDa
Shipping	Ice bag