


Anti-Eif3B Polyclonal Antibody For Wb, If/Icc, Elisa - P55884

Item Number: CM0017776



Introduction

Rabbit polyclonal antibody against human eIF3B, a subunit of the translation initiation factor eIF-3 complex. Suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Predicted molecular weight ~93 kDa. Ideal for translation research.

[Learn More](#)

Parameter	Value
Product Name	eIF3B Rabbit pAb
Remarks/Alias	PRT1, EIF3S9, EIF3-ETA, EIF3-P110, EIF3-P116, eIF3B
Species	Human
GeneID (Human)	8662
GeneID	8662
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-350 of human eIF3B (NP_003742.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55884
Protein Weight	93kDa
Shipping	Ice bag