

Anti-Eif3B Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P55884

Item Number: CM0012480



Introduction

Rabbit monoclonal antibody against human eIF3B, validated for Western blot, immunofluorescence (IF/ICC), and ELISA. Cross-reacts with mouse and rat. Suitable for translation initiation and cell proliferation research.

[Learn More](#)

Parameter	Value
Product Name	eIF3B Rabbit mAb
Remarks/Alias	PRT1; EIF3S9; EIF3-ETA; EIF3-P110; EIF3-P116; eIF3B
Species	Human
GeneID (Human)	8662
GeneID	8662
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 101-200 of human eIF3B (NP_003742.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55884
Protein Weight	93kDa
Shipping	Ice bag