


Anti-Eif3C Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q99613

Item Number: CM0030991



Introduction

EIF3C Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat EIF3C (eIF3-p110), a core subunit of the eIF-3 translation initiation complex. Ideal for protein synthesis and cell proliferation studies.

[Learn More](#)

Parameter	Value
Product Name	EIF3C Rabbit pAb
Remarks/Alias	EIF3CL; EIF3S8; eIF3-p110; EIF3C
Species	Human
Gene ID (Human)	8663
Gene ID	8663
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 694-913 of human EIF3C (NP_003743.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99613
Protein Weight	106kDa
Shipping	Ice bag