


Anti-Eif3D Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O15371

Item Number: CM0025185



Introduction

Rabbit polyclonal antibody against human EIF3D (eIF3-zeta). Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Suitable for studying translation initiation, cell proliferation, and host-virus interactions.

[Learn More](#)

Parameter	Value
Product Name	EIF3D Rabbit pAb
Remarks/Alias	EIF3S7; eIF3-p66; eIF3-zeta; EIF3D
Species	Human
Gene ID (Human)	8664
Gene ID	8664
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15371
Protein Weight	64kDa
Shipping	Ice bag