


Anti-Eif3E Monoclonal Antibody For Wb, Ihc-P, Elisa - P60228

Item Number: CM0009258



Introduction

Rabbit monoclonal antibody against eIF3e (EIF3S6/INT6) for Western blot, IHC-P, and ELISA. Validated in human, mouse, and rat. Suitable for studying translation initiation and NMD.

[Learn More](#)

Parameter	Value
Product Name	eIF3e Rabbit mAb
Remarks/Alias	INT6, EIF3S6, EIF3-P48, eIF3-p46, eIF3e
Species	Human
GeneID (Human)	3646
GeneID	3646
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-350 of human eIF3e (P60228)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60228
Protein Weight	52kDa
Shipping	Ice bag