


Anti-Eif3E Rabbit Polyclonal Antibody For Wb, Elisa - P60228

Item Number: CM0026553



Introduction

CamelBio offers a high-quality rabbit polyclonal antibody targeting eIF3e (INT6, EIF3S6). Validated for Western blot and ELISA in human and mouse samples, this antibody supports research on protein translation initiation and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	eIF3e Rabbit pAb
Remarks/Alias	INT6; EIF3S6; EIF3-P48; eIF3-p46; eIF3e
Species	Human
GeneID (Human)	3646
GeneID	3646
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-445 of human eIF3e (NP_001559.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P60228
Protein Weight	52kDa
Shipping	Ice bag