


Anti-Eif3G Polyclonal Antibody For Wb, If/Icc, Elisa - 075821

Item Number: CM0030496



Introduction

High-quality rabbit polyclonal antibody targeting human EIF3G, validated for Western blot, immunofluorescence/immunocytochemistry, and ELISA. EIF3G is a core component of the eukaryotic translation initiation factor 3 complex, essential for protein synthesis and cell cycle regulation. This antibody cross-reacts with mouse and rat samples, providing versatile utility in translational research.

[Learn More](#)

Parameter	Value
Product Name	EIF3G Rabbit pAb
Remarks/Alias	EIF3S4; EIF3-P42; eIF3-p44; eIF3-delta; EIF3G
Species	Human
Gene ID (Human)	8666
Gene ID	8666
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of human EIF3G (NP_003746.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075821
Protein Weight	36kDa
Shipping	Ice bag