


Anti-Eif3J Polyclonal Antibody For Wb, If/Icc, Elisa - 075822

Item Number: CM0030390



Introduction

EIF3J Rabbit pAb (075822) is a polyclonal antibody validated for WB, IF/ICC, and ELISA in human, mouse, and rat. Ideal for research on translation initiation and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	EIF3J Rabbit pAb
Remarks/Alias	EIF3S1; eIF3-p35; eIF3-alpha; EIF3J
Species	Human
Gene ID (Human)	8669
Gene ID	8669
Immunogen	Recombinant protein corresponding to amino acids 1-258 of human EIF3J (NP_003749.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075822
Protein Weight	29kDa
Shipping	Ice bag