


Anti-Eif3K Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Ubq5

Item Number: CM0025506



Introduction

High-quality rabbit polyclonal antibody against EIF3K, validated for Western blot, IHC-P, and ELISA in human, mouse, and rat. Ideal for studying translation initiation and protein synthesis.

[Learn More](#)

Parameter	Value
Product Name	EIF3K Rabbit pAb
Remarks/Alias	M9; ARG134; PLAC24; PTD001; EIF3S12; HSPC029; MSTP001; PLAC-24; PRO1474; EIF3-p28; EIF3K
Species	Human
GeneID (Human)	27335
GeneID	27335
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-218 of human EIF3K (NP_037366.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBQ5
Protein Weight	25kDa
Shipping	Ice bag