


Anti-Eif3L Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y262

Item Number: CM0028040



Introduction

Anti-EIF3L Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Component of eIF-3 translation initiation complex, involved in protein synthesis and cell proliferation. Ideal for IVD and research use.

[Learn More](#)

Parameter	Value
Product Name	EIF3L Rabbit pAb
Remarks/Alias	EIF3EIP; EIF3S11; HSPC021; HSPC025; MSTP005; EIF3S6IP; EIF3L
Species	Human
GeneID (Human)	51386
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human EIF3L (NP_057175.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y262
Protein Weight	67kDa
Shipping	Ice bag