


Anti-Eif3M Polyclonal Antibody For Wb, Ihc-P, Elisa - Q7L2H7

Item Number: CM0028126



Introduction

Rabbit polyclonal antibody against human EIF3M (eIF3m, hFL-B5), validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. EIF3M is a subunit of the eIF-3 complex essential for translation initiation, relevant in cancer and viral infection research.

[Learn More](#)

Parameter	Value
Product Name	EIF3M Rabbit pAb
Remarks/Alias	B5; GA17; PCID1; TANGO7; hfl-B5; EIF3M
Species	Human
Gene ID (Human)	10480
Gene ID	10480
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-374 of human EIF3M (NP_006351.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7L2H7
Protein Weight	42kDa
Shipping	Ice bag