

Anti-Eif4A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P60842 / Q14240 / P38919

Item Number: CM0020964



Introduction

eIF4A rabbit polyclonal antibody for WB, IF/ICC, ELISA. Recognizes human, mouse, rat eIF4A (P60842/Q14240/P38919). Ideal for translation initiation and RNA helicase research.

[Learn More](#)

Parameter	Value
Product Name	eIF4A Rabbit pAb
Species	Human
Gene ID	1973/1974/9775
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human eIF4A (NP_001407.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60842/Q14240/P38919
Shipping	Ice bag