

Anti-Eif4A1 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P60842

Item Number: CM0028795



Introduction

eIF4A1 Rabbit pAb validated for WB, IF/ICC, IP, ELISA; reacts with human, mouse, rat. Ideal for studying translation initiation and RNA helicase function in research and diagnostics.

[Learn More](#)

Parameter	Value
Product Name	eIF4A1 Rabbit pAb
Remarks/Alias	DDX2A; EIF4A; EIF-4A; eIF4A-I; eIF-4A-I; eIF4A1
Species	Human
Gene ID	1973 (human)
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human eIF4A1 (NP_001407.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60842
Protein Weight	46kDa
Shipping	Ice bag