

Anti-Eif4A1/Eif4A2/Eif4A3 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P60842/Q14240/P38919

Item Number: CM0007538



Introduction

Recombinant rabbit monoclonal antibody targeting human EIF4A1/EIF4A2/EIF4A3 with cross-reactivity to mouse and rat. Validated for WB, IHC-P, IF/ICC, IP, and ELISA. Ideal for studying translation initiation and RNA helicase function.

[Learn More](#)

Parameter	Value
Product Name	EIF4A1/EIF4A2/EIF4A3 Rabbit mAb
Remarks/Alias	DDX2A; EIF4A; EIF-4A; eIF4A-I; eIF-4A-I; DDX2B; EIF4F; BM-010; eIF4A-II; e-IF-4A-II; Fal1; RCPS; DDX48; MUK34; NUK34; NMP265; eIF4AIII; eIF4A-III; e-IF-4A-III; EIF4A1/EIF4A2/EIF4A3
Species	Human
GeneID (Human)	1973/1974/9775
GeneID	1973/1974/9775
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human EIF4A1, EIF4A2, and EIF4A3
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60842/Q14240/P38919
Protein Weight	46kDa
Shipping	Ice bag