


Anti-Eif4A2 Polyclonal Antibody For Wb, Elisa - Q14240

Item Number: CM0025036



Introduction

High-quality anti-EIF4A2 rabbit polyclonal antibody validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Suitable for studying translation initiation and RNA helicase activity.

[Learn More](#)

Parameter	Value
Product Name	EIF4A2 Rabbit pAb
Remarks/Alias	DDX2B; EIF4A; EIF4F; BM-010; eIF4A-II; eIF-4A-II; EIF4A2
Species	Human
GeneID (Human)	1974
GeneID	1974
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human EIF4A2 (NP_001958.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14240
Protein Weight	46kDa
Shipping	Ice bag