

Anti-Eif4A3 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P38919

Item Number: CM0005751



Introduction

Anti-eIF4A3 rabbit monoclonal antibody validated for WB, IF/ICC, and ELISA. Recognizes human, mouse, and rat eIF4A3 protein (P38919). Ideal for studying splicing and EJC dynamics.

[Learn More](#)

Parameter	Value
Product Name	eIF4A3 Rabbit mAb
Remarks/Alias	Fal1; RCPS; DDX48; MUK34; NUK34; NMP265; eIF4AIII; eIF4A-III; eIF-4A-III; eIF4A3
Species	Human
Gene ID	9775 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-125 of human eIF4A3 (P38919)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P38919
Protein Weight	47kDa
Shipping	Ice bag