

Anti-Eif4A3 Rabbit Polyclonal Antibody For Wb, Elisa - P38919

Item Number: CM0025711



Introduction

Rabbit polyclonal antibody targeting human eIF4A3 (DDX48/NMP265), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for studying mRNA splicing, EJC assembly, and nonsense-mediated decay.

[Learn More](#)

Parameter	Value
Product Name	eIF4A3 Rabbit pAb
Remarks/Alias	Fal1; RCPS; DDX48; MUK34; NUK34; NMP265; eIF4AIII; eIF4A-III; eIF-4A-III; eIF4A3
Species	Human
GeneID (Human)	9775
GeneID	9775
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human eIF4A3 (NP_055555.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P38919
Protein Weight	47kDa
Shipping	Ice bag