


Anti-Eif4B Polyclonal Antibody For Wb, If/Icc, Elisa - P23588

Item Number: CM0017067



Introduction

Rabbit polyclonal antibody against human EIF4B, validated for WB, IF/ICC, and ELISA. Detects human, mouse, rat EIF4B; suitable for translation initiation studies.

[Learn More](#)

Parameter	Value
Product Name	EIF4B Rabbit pAb
Remarks/Alias	EIF-4B; PRO1843; EIF4B
Species	Human
GeneID (Human)	1975
GeneID	1975
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 352-611 of human EIF4B (NP_001408.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23588
Protein Weight	69kDa
Shipping	Ice bag