


Anti-Eif4B Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P23588

Item Number: CM0026672



Introduction

EIF4B Rabbit pAb is a polyclonal antibody targeting human and mouse eukaryotic translation initiation factor 4B. Validated for WB, IF/ICC, IP, and ELISA, it supports research on protein biosynthesis and translational control mechanisms.

[Learn More](#)

Parameter	Value
Product Name	EIF4B Rabbit pAb
Remarks/Alias	EIF-4B; PRO1843; EIF4B
Species	Human
Gene ID (Human)	1975
Gene ID	1975
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 352-611 of human EIF4B (NP_001408.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P23588
Protein Weight	69kDa
Shipping	Ice bag