

Anti-Eif4E Polyclonal Antibody For Wb, If/Icc, Elisa - P06730

Item Number: CM0013290



Introduction

Rabbit polyclonal antibody against human eIF4E, validated for WB, IF/ICC, ELISA. Recognizes the 25 kDa translation initiation factor involved in mRNA cap binding, nuclear export, and translational control.

[Learn More](#)

Parameter	Value
Product Name	eIF4E Rabbit pAb
Remarks/Alias	CBP; EIF4F; AUTS19; EIF4E1; eIF-4E; EIF4EL1; eIF4E
Species	Human
Gene ID (Human)	1977
Gene ID	1977
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-217 of human eIF4E (NP_001959.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P06730
Protein Weight	25kDa
Shipping	Ice bag