

## Anti-Eif4E Rabbit Monoclonal Antibody [Kd Validated] For Wb, If/Icc, Ip, Elisa - P06730

Item Number: CM0009466



### Introduction

High-quality rabbit monoclonal antibody against human eIF4E, validated for WB, IF/ICC, IP, ELISA. Cross-reacts with mouse, rat, monkey. Ideal for translation initiation and mRNA export studies - P06730.

[Learn More](#)

| Parameter        | Value                                                                                               |
|------------------|-----------------------------------------------------------------------------------------------------|
| Product Name     | [KD Validated] eIF4E Rabbit mAb                                                                     |
| Remarks/Alias    | CBP; EIF4F; AUTS19; EIF4E1; eIF-4E; EIF4EL1                                                         |
| Species          | Human                                                                                               |
| Gene ID          | 1977                                                                                                |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 1-100 of human eIF4E (NP_001959.1) |
| Source           | Rabbit                                                                                              |
| Category         | Monoclonal Antibodies                                                                               |
| Application      | WB, IF/ICC, IP, ELISA                                                                               |
| Cross Reactivity | Human, Mouse, Rat, Monkey                                                                           |
| SWISS            | P06730                                                                                              |
| Protein Weight   | 25kDa/27kDa/29kDa                                                                                   |
| Shipping         | Ice bag                                                                                             |