

Anti-Eif4E2 Polyclonal Antibody For Wb - O60573

Item Number: CM0024764



Introduction

Rabbit polyclonal antibody targeting human EIF4E2 (4EHP), validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes Human, Mouse, Rat. Key translational repressor involved in miRNA-mediated silencing and viral response.

[Learn More](#)

Parameter	Value
Product Name	EIF4E2 Rabbit pAb
Remarks/Alias	4EHP, IF4e, 4E-LP, h4EHP, EIF4EL3, EIF4E2
Species	Human
GeneID (Human)	9470
GeneID	9470
Immunogen	Recombinant fusion protein containing amino acids 1-240 of human EIF4E2 (NP_004837.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60573
Protein Weight	28kDa
Shipping	Ice bag