


Anti-Eif4Ebp1 (4E-Bp1) Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13541

Item Number: CM0004293



Introduction

Rabbit monoclonal antibody against eIF4EBP1 (4E-BP1, PHAS-I) validated for WB, IHC-P, IF/ICC, ELISA. Suitable for human, mouse, rat studies of translation regulation, mTOR signaling, and cancer research. Uniprot: Q13541.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] eIF4EBP1 Rabbit mAb
Remarks/Alias	BP-1; 4EBP1; 4E-BP1; PHAS-I; eIF4EBP1
Species	Human
Gene ID (Human)	1978
Gene ID	1978
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 20-118 of human eIF4EBP1(NP_004086.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13541
Protein Weight	13kDa
Shipping	Ice bag