

## Anti-Eif4Ebp1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13541

Item Number: CM0013131



### Introduction

Rabbit polyclonal antibody against human eIF4EBP1, validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples. Ideal for translation regulation studies.

[Learn More](#)

| Parameter        | Value                                                                                                                         |
|------------------|-------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | eIF4EBP1 Rabbit pAb                                                                                                           |
| Remarks/Alias    | BP-1; 4EBP1; 4E-BP1; PHAS-I; eIF4EBP1                                                                                         |
| Species          | Human                                                                                                                         |
| Gene ID          | 1978                                                                                                                          |
| Immunogen        | Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 50-118 of human eIF4EBP1 (NP_004086.1). |
| Source           | Rabbit                                                                                                                        |
| Category         | Polyclonal Antibodies                                                                                                         |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                                                                             |
| SWISS            | Q13541                                                                                                                        |
| Protein Weight   | 13kDa                                                                                                                         |
| Shipping         | Ice bag                                                                                                                       |