

Anti-Eif4Ebp1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q13541

Item Number: CM0007479



Introduction

High specificity anti-eIF4EBP1 rabbit monoclonal antibody validated for WB, IF/ICC, and ELISA; detects human eIF4EBP1 (13 kDa) for translation regulation studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] eIF4EBP1 Rabbit mAb
Remarks/Alias	BP-1; 4EBP1; 4E-BP1; PHAS-I; [KD Validated] eIF4EBP1
Species	Human
GeneID (Human)	1978
GeneID	1978
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q13541
Protein Weight	13kDa
Shipping	Ice bag