


Anti-Eif4Ebp1 Rabbit Polyclonal Antibody For Ihc-P And Elisa - Q13541

Item Number: CM0019761



Introduction

Rabbit polyclonal antibody against human eIF4EBP1 (4EBP1) for IHC-P and ELISA, cross-reacts with mouse and rat. Ideal for investigating mTORC1 pathway and translation initiation control.

[Learn More](#)

Parameter	Value
Product Name	eIF4EBP1 Rabbit pAb
Remarks/Alias	BP-1; 4EBP1; 4E-BP1; PHAS-I; eIF4EBP1
Species	Human
GeneID (human)	1978
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-118 of human eIF4EBP1 (NP_004086.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13541
Protein Weight	13kDa
Shipping	Ice bag