


Anti-Eif4G1 Polyclonal Antibody For Wb And Elisa - Q04637

Item Number: CM0017349



Introduction

Rabbit polyclonal antibody against human EIF4G1 (eIF4G, p220), validated for WB and ELISA. Suitable for studying translation initiation, stress granules, and Parkinson's disease-related pathways.

[Learn More](#)

Parameter	Value
Product Name	EIF4G1 Rabbit pAb
Remarks/Alias	P220, EIF4F, EIF4G, EIF4GI, PARK18, EIF-4G1, EIF4G1
Species	Human
GeneID (Human)	1981
GeneID	1981
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 550-650 of human EIF4G1 (NP_886553.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q04637
Protein Weight	175kDa
Shipping	Ice bag