


Anti-Eif4G1 Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q04637

Item Number: CM0027562



Introduction

Rabbit polyclonal antibody against human EIF4G1 (eIF4G, p220), validated for WB, IHC-P, IP, ELISA. Ideal for studying translation initiation and stress response pathways.

[Learn More](#)

Parameter	Value
Product Name	EIF4G1 Rabbit pAb
Remarks/Alias	P220; EIF4F; EIF4G; EIF4GI; PARK18; EIF-4G1; EIF4G1
Species	Human
Gene ID (Human)	1981
Gene ID	1981
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 400-500 of human EIF4G1 (NP_886553.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human
SWISS	Q04637
Protein Weight	175kDa
Shipping	Ice bag