

Anti-Eif4G2/P97 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P78344

Item Number: CM0028671



Introduction

Rabbit polyclonal antibody to EIF4G2/p97, validated for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Suitable for translation initiation research.

[Learn More](#)

Specification	Details
Product Name	EIF4G2/p97 Rabbit pAb
Remarks/Alias	P97; AAG1; DAP5; NAT1; EIF4G2/p97
Species	Human
Gene ID (Human)	1982
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 320-490 of human EIF4G2/p97 (NP_001409.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78344
Protein Weight	102kDa
Shipping	Ice bag