


Anti-Eif5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P55010

Item Number: CM0025280



Introduction

This rabbit polyclonal antibody targets EIF5, suitable for WB, IF/ICC, and ELISA in human, mouse, and rat samples. EIF5 is a GTPase-activating protein critical for 43S pre-initiation complex function in mRNA translation initiation.

[Learn More](#)

Parameter	Value
Product Name	EIF5 Rabbit pAb
Remarks/Alias	EIF-5; EIF-5A; EIF5
Species	Human
GeneID	1983
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-235 of human EIF5 (NP_892116.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55010
Protein Weight	49kDa
Shipping	Ice bag