

Anti-Eif5A Rabbit Monoclonal Antibody For Wb, Elisa - P63241

Item Number: CM0010790



Introduction

High-quality anti-eIF5A rabbit monoclonal antibody validated for Western blot (WB) and ELISA applications. Cross-reacts with human, mouse, and rat. Ideal for translation elongation, ribosome stalling, and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	eIF5A Rabbit mAb
Remarks/Alias	FABAS; EIF-5A; EIF5A1; eIF-4D; eIF5AI; eIF5A
Species	Human
GeneID (human)	1984
GeneID	1984
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human eIF5A (P63241).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63241
Protein Weight	17kDa
Shipping	Ice bag