

Anti-Eif5A2 Polyclonal Antibody For Wb, Elisa - Q9Gzv4

Item Number: CM0026354

Introduction

Rabbit polyclonal antibody against human EIF5A2, validated for WB and ELISA. Detects human, mouse, rat EIF5A2. Ideal for translation and cancer research. SWISS: Q9GZV4.

[Learn More](#)

Parameter	Value
Product Name	EIF5A2 Rabbit pAb
Remarks/Alias	EIF-5A2; eIF5AII; EIF5A2
Species	Human
GeneID (Human)	56648
GeneID	56648
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-153 of human EIF5A2 (NP_065123.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZV4
Protein Weight	17kDa
Shipping	Ice bag