


## Anti-Eif5B Polyclonal Antibody For Wb, If/Icc, Elisa - O60841

Item Number: CM0014033



### Introduction

Rabbit polyclonal antibody targeting human EIF5B (IF2) validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for translation initiation studies and ribosomal subunit joining research.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | EIF5B Rabbit pAb                                                                                                                      |
| Remarks/Alias    | IF2; EIF5B                                                                                                                            |
| Species          | Human                                                                                                                                 |
| GeneID (human)   | 9669                                                                                                                                  |
| GeneID           | 9669                                                                                                                                  |
| Immunogen        | Recombinant Protein, Recombinant fusion protein containing a sequence corresponding to amino acids 1-273 of human EIF5B (NP_056988.3) |
| Source           | Rabbit                                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                                 |
| Application      | WB, IF/ICC, ELISA                                                                                                                     |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                     |
| SWISS            | O60841                                                                                                                                |
| Protein Weight   | 139kDa                                                                                                                                |
| Shipping         | Ice bag                                                                                                                               |