

Anti-Elastin (Eln) Rabbit Monoclonal Antibody For Wb, Ihc-P, If-P, Elisa - P15502

Item Number: CM0007453

Introduction

Rabbit monoclonal antibody against human elastin (ELN), validated for WB, IHC-P, IF-P, and ELISA. Cross-reacts with human and rat. Targets amino acids 27-200. For extracellular matrix studies.

[Learn More](#)

Parameter	Value
Product Name	ELN Rabbit mAb
Remarks/Alias	WS; WBS; SVAS; ADCL1; ELN
Species	Human
Gene ID	2006 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-200 of human ELN (NP_001265844.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P15502
Protein Weight	68kDa
Shipping	Ice bag