

Anti-Elastin Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P15502

Item Number: CM0021316

Introduction

Anti-Elastin (ELN) rabbit polyclonal antibody validated for Western blot, immunofluorescence, immunohistochemistry, and ELISA; reacts with human, mouse, rat; suitable for extracellular matrix and vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	ELN Rabbit pAb
Remarks/Alias	WS; WBS; SVAS; ADCL1; ELN
Species	Human
GeneID	2006
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-200 of human ELN (NP_001265868.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15502
Protein Weight	68 kDa
Shipping	Ice bag