

Anti-Elavl2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q12926

Item Number: CM0030729

Introduction

Rabbit polyclonal antibody against human ELAVL2 (HuB, HEL-N1), validated in WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Suitable for neuroscience and RNA-binding protein research.

[Learn More](#)

Parameter	Value
Product Name	ELAVL2 Rabbit pAb
Remarks/Alias	HUB; HELN1; HEL-N1; ELAVL2
Species	Human
GeneID (Human)	1993
GeneID	1993
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human ELAVL2 (NP_004423.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12926
Protein Weight	40kDa
Shipping	Ice bag