

Anti-Elavl3 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q14576

Item Number: CM0028986

Introduction

Rabbit polyclonal antibody targeting human ELAVL3 (HuC, PLE21). Validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for neuroscience research on neuronal differentiation and mRNA regulation. - 40kDa, Q14576.

[Learn More](#)

Parameter	Value
Product Name	ELAVL3 Rabbit pAb
Remarks/Alias	HUC; HUCL; PLE21; ELAVL3
Species	Human
Gene ID (Human)	1995
Gene ID	1995
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 161-260 of human ELAVL3 (NP_001411.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14576
Protein Weight	40kDa
Shipping	Ice bag