

Anti-Elavl4 Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P26378

Item Number: CM0016464



Introduction

Rabbit polyclonal antibody targeting human ELAVL4 (HuD) for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for neurobiology and paraneoplastic syndrome research. Predicted MW 42 kDa. SWISS P26378.

[Learn More](#)

Parameter	Value
Product Name	ELAVL4 Rabbit pAb
Remarks/Alias	HUD; PNEM; ELAVL4
Species	Human
GeneID	1996
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human ELAVL4 (NP_068771.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P26378
Protein Weight	42kDa
Shipping	Ice bag