

Anti-Elf2 Rabbit Polyclonal Antibody For Wb And Elisa - Q15723

Item Number: CM0025470

Introduction

High-quality rabbit polyclonal antibody against human ELF2 (EU32/NERF). Validated for WB and ELISA. Cross-reacts with human and mouse ELF2. Target protein weight 64 kDa.

[Learn More](#)

Parameter	Value
Product Name	ELF2 Rabbit pAb
Remarks/Alias	EU32; NERF; NERF-2; NERF-1A; NERF-1B; NERF-1a; b; ELF2
Species	Human
Gene ID (Human)	1998
Gene ID	1998
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15723
Protein Weight	64kDa
Shipping	Ice bag