

Anti-Elf4 Polyclonal Antibody For If/Icc And Elisa - Q99607

Item Number: CM0025948

Introduction

Rabbit polyclonal antibody against human ELF4, validated for immunofluorescence (IF/ICC) and ELISA. Recombinant immunogen, cross-reacts with human ELF4 (UniProt Q99607). Suitable for studying transcriptional regulation, innate immunity, and inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	ELF4 Rabbit pAb
Remarks/Alias	MEF; ELFR; AIFBL2; ELF4
Species	Human
Gene ID (Human)	2000
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human ELF4 (NP_001412.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q99607
Protein Weight	71kDa
Shipping	Ice bag