


Anti-Elf5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Ukw6

Item Number: CM0025159



Introduction

Rabbit polyclonal antibody against ELF5, validated in WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Suitable for epithelial differentiation and transcription regulation studies. UniProt Q9UKW6.

[Learn More](#)

Parameter	Value
Product Name	ELF5 Rabbit pAb
Remarks/Alias	ESE2; ELF5
Species	Human
Gene ID (Human)	2001
Gene ID	2001
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-150 of human ELF5 (NP_938195.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKW6
Protein Weight	31kDa
Shipping	Ice bag