

Anti-Elk1 Polyclonal Antibody For Wb, Elisa - P19419

Item Number: CM0016139



Introduction

ELK1 Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with human, mouse, rat. Targets the C-terminal region of the ETS-domain transcription factor ELK1, involved in MAPK signaling and immediate early gene induction.

[Learn More](#)

Parameter	Value
Product Name	ELK1 Rabbit pAb
Remarks/Alias	ELK1
Species	Human
Gene ID	2002
Immunogen	Recombinant fusion protein containing amino acids 201-297 of human ELK1 (NP_005220.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19419
Protein Weight	45kDa
Shipping	Ice bag