

Anti-Elk1 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - P19419

Item Number: CM0003497



Introduction

Rabbit monoclonal ELK1 antibody validated for WB, IP, and ELISA in human, mouse, and rat. Targets the ETS domain transcription factor Elk-1, a key mediator of MAPK signaling and immediate-early gene activation. Research use only.

[Learn More](#)

Parameter	Value
Product Name	ELK1 Rabbit mAb
Remarks/Alias	ELK1
Species	Human
Gene ID (Human)	2002
Gene ID	2002
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 329-428 of human ELK1 (NP_001107595.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19419
Protein Weight	45kDa
Shipping	Ice bag