

Anti-ELL Rabbit Polyclonal Antibody For Wb, Elisa - P55199

Item Number: CM0023643



Introduction

Rabbit polyclonal antibody against human ELL (Eleven-nineteen lysine-rich leukemia protein), validated for Western blot and ELISA, cross-reacts with mouse. Suitable for studying RNA polymerase II elongation and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	ELL Rabbit pAb
Remarks/Alias	MEN; ELL1; PPP1R68; C19orf17; ELL
Species	Human
GeneID	8178
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 350-450 of human ELL (NP_006523.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P55199
Protein Weight	68kDa
Shipping	Ice bag