


Anti-ELL2 Rabbit Polyclonal Antibody For Wb And Elisa - 000472

Item Number: CM0024474



Introduction

High-quality anti-ELL2 rabbit polyclonal antibody validated for WB and ELISA. Cross-reacts with rat. Ideal for investigating RNA polymerase II elongation, snRNA transcription, and immunoglobulin processing.

[Learn More](#)

Parameter	Value
Product Name	ELL2 Rabbit pAb
Remarks/Alias	MRCCAT1; ELL2
Species	Human
GeneID (Human)	22936
GeneID	22936
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-292 of human ELL2 (NP_036213.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	000472
Protein Weight	72kDa
Shipping	Ice bag