

Anti-Elmo1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q92556

Item Number: CM0005685



Introduction

High-quality rabbit monoclonal antibody targeting human ELMO1 (Q92556) for WB, IHC-P, ELISA. Validated across human, mouse, and rat. Ideal for studying cytoskeletal rearrangements, phagocytosis, and cell motility.

[Learn More](#)

Parameter	Value
Product Name	ELMO1 Rabbit mAb
Remarks/Alias	CED12; CED-12; ELMO-1; ELMO1
Species	Human
GeneID (human)	9844
GeneID	9844
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 628-727 of human ELMO1 (Q92556).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92556
Protein Weight	84kDa
Shipping	Ice bag