


Anti-Elmo1 Rabbit Polyclonal Antibody For Wb, Elisa - Q92556

Item Number: CM0013321



Introduction

High-quality rabbit polyclonal antibody targeting human ELMO1, validated for WB and ELISA. Reacts with human and mouse samples. Ideal for studying phagocytosis, cell motility, and DOCK-GEF/RAC1 signaling.

[Learn More](#)

Parameter	Value
Product Name	ELMO1 Rabbit pAb
Remarks/Alias	CED12; CED-12; ELMO-1; ELMO1
Species	Human
GeneID (Human)	9844
GeneID	9844
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 527-727 of human ELMO1 (NP_055615.8)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92556
Protein Weight	84kDa
Shipping	Ice bag