

Anti-Elmo2 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Jj3

Item Number: CM0015538



Introduction

Rabbit polyclonal antibody targeting human ELMO2, validated for WB and ELISA. Cross-reacts with human, mouse, rat. Useful for phagocytosis, cell motility, and engulfment studies. Immunogen: recombinant fusion protein (aa 42-120).

[Learn More](#)

Parameter	Value
Product Name	ELMO2 Rabbit pAb
Remarks/Alias	VMPI; CED12; CED-12; ELMO-2; Ced-12A; ELMO2
Species	Human
Gene ID	63916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-120 of human ELMO2 (NP_877496.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96Jj3
Protein Weight	83kDa
Shipping	Ice bag