

Anti-Elmo3 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Bj8

Item Number: CM0025641

Introduction

Human ELMO3 rabbit polyclonal antibody validated for WB and ELISA; cross-reacts with human, mouse, rat. Target is a phagocytosis and cell motility protein involved in cytoskeletal rearrangements. Ideal for apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	ELMO3 Rabbit pAb
Remarks/Alias	CED12; CED-12; ELMO-3; ELMO3
Species	Human
GeneID (Human)	79767
GeneID	79767
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 504-773 of human ELMO3 (NP_078988.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96Bj8
Protein Weight	81kDa
Shipping	Ice bag